# Information for Patients

## Anatomical Breast Implants

|                      |                      |
|----------------------|----------------------|
| <b>EN</b> English    | <b>RU</b> Russian    |
| <b>FR</b> French     | <b>CS</b> Czech      |
| <b>DE</b> German     | <b>ET</b> Estonian   |
| <b>IT</b> Italian    | <b>HU</b> Hungarian  |
| <b>ES</b> Spanish    | <b>PL</b> Polish     |
| <b>PT</b> Portuguese | <b>SL</b> Slovenian  |
| <b>DA</b> Danish     | <b>SV</b> Swedish    |
| <b>NL</b> Dutch      | <b>LT</b> Lithuanian |
| <b>EL</b> Greek      | <b>BG</b> Bulgarian  |
| <b>LV</b> Latvian    | <b>HR</b> Croatian   |
| <b>SK</b> Slovak     | <b>FI</b> Finnish    |
| <b>RO</b> Romanian   | <b>NO</b> Norwegian  |
| <b>TR</b> Turkish    | <b>IS</b> Icelandic  |

| A Confident Choice for Life™

## EN Information for patients

### Implant Card Purpose

The implant card filled by the health care professional should be provided to you after the surgery.

The implant card must be kept on your person at all times and enables you:

- To identify the implanted devices and to get access to other information related to the implanted device (e.g. via EUDAMED, and other websites).
- To identify yourselves as persons requiring special care in relevant situations e.g. medical examinations, security checks.
- To inform emergency clinical staff or first responder about special care/needs for you in case of emergency situations.

### The International Implant Card is linked with the GCA Comfort™ Plus Warranty.

This card applies to all GC Aesthetics® Silicone Gel Filled Breast Implants (certified under MDR 2017/745) implanted by an appropriately qualified and licensed Physician in accordance with accepted surgical procedures and techniques, and in strict accordance with the GCA® "Instruction for Use" supplied with the product.

### Patient leaflet Purpose

The purpose of this leaflet is to provide additional information you should be aware in terms of:

1. Information about the **expected lifetime of the device** and any necessary follow up.
2. Any information to **ensure safe use of the device** by the patient, including information on the materials and substances to which you can be exposed.
3. **Warnings, precautions or measures** to be taken by the patient or a healthcare professional with regard to interference with any possible external influences, medical examinations or environmental conditions.

The updated version of this leaflet is available on our information website for patients, by the link detailed on the implant card:  
<https://www.gcaesthetics.com/patient>

### 1. Product Lifetime

Breast implants are not lifetime devices. Be aware that they may need to be removed or replaced at some future date which can lead to revision surgery. Many surgeons suggest replacement after 10-20 years. But unless there is a problem, this may not be necessary.

### 2. Necessary Follow-Up / Information to Ensure Safety

#### Postoperative Care

Your surgeon will provide advice on how to take care of yourself after you have breast implant surgery. These may include you wearing a supportive bra 24 hours a day or changing the adhesive bandages when necessary. Your surgeon may also prescribe prophylactic antibiotics. These and other measures are at your surgeons opinion and are based on their personal experience and overall treatment plan.

Not following the surgeon's recommendations may lead to complications.

You should consult your surgeon if you suspect any complications. In the case of trauma to the area, please speak to your surgeon. Trauma can include extreme massaging of the breast regions, by some sporting activity or by using seat belts. Contact sports are not recommended. Tanning at salons or sunbathing can influence scar healing. Direct scar exposure to the sun should be avoided for about six months after surgery. It is recommended that you consult a physician or pharmacist before the use of topical medicines (e.g., steroids) in the breast area.

If any surgery of the breast area is scheduled, you should inform the physician or the surgeon of the presence of an implant and carry your implant card.

The postoperative care recommendations presented within this section of the patient brochure does not include all possible options. Your surgeon should provide postoperative care details.

#### After Surgery

Once breast augmentation surgery has taken place, a follow-up appointment with a surgeon will be scheduled. This is routine and will be to check on your general health after the operation.

Following this, contact with a surgeon should be limited if you continue to feel fit and healthy.

However, you should always remember to consult a surgeon, physician or pharmacist before the use of topical medicines (e.g. steroids) in the breast area. You should also inform them of the presence of implants if you require any surgery in the future.

#### Self-Examination and Auto palpation

It's recommended that adult women of all ages perform a breast self-exam at least once a month. A regular self-exam will help you to be familiar with how your breasts look and feel. You can contact a healthcare professional if you notice any changes.

One of the ways you can perform a breast self-exam is in front of a mirror by looking for any changes in contour, any swelling, or dimpling of the skin, or changes in the nipples. You should do this with your arms at your sides and your arms high overhead. Next, rest your palms on your hips and press firmly to flex your chest muscles. Look for any dimpling, puckering, or changes, particularly on one side.

#### Breast and Nipple Piercing Procedures

Women with breast implants who undergo body piercing procedures to the breast region must consider the possibility that an infection could develop following this procedure. Treatment including antibiotics, possible removal of the implant, or additional surgery may be necessary.

#### Materials

These medical devices are manufactured from medical grade silicone materials which are suitable for long term implantation. Master Files have been filed with the U.S FDA.

The following three tables provide quantitative and qualitative information on the materials and substances which patients can be exposed based on chemical characterization of representative devices. The materials and device have been subjected to biocompatibility testing and evaluation and risk assessments to demonstrate their biological safety. However, individual responses to chemicals may vary, and all reactions cannot be predicted.

Table 1. Breast Implant Device Materials

| Device materials     | Implant component         |
|----------------------|---------------------------|
| Silicone dispersions | Shell                     |
| Silicone adhesive    | Shell                     |
| Silicone gel         | Gel fill                  |
| Silicone dispersions | Closure Patch             |
| Silicone Elastomer   | Orientation markers       |
| Platinum catalyst    | Shell and Gel fill        |
| Organotin catalyst   | Shell (silicone adhesive) |

Table 2. Chemicals Released by Breast Implants

| Compounds                                                                     | Whole device (µg/g or ppm) | Compounds                           | Whole device (µg/g or ppm) |
|-------------------------------------------------------------------------------|----------------------------|-------------------------------------|----------------------------|
| <b>Volatiles</b> - Chemicals that can be released by breast implants as a gas |                            |                                     |                            |
| Benzene                                                                       | Not Detected               | Ethylbenzene                        | Not Detected               |
| Bromobenzene                                                                  | Not Detected               | Hexachlorobutadiene                 | Not Detected               |
| Bromodichloromethane                                                          | Not Detected               | Isopropylbenzene                    | Not Detected               |
| Bromoform                                                                     | Not Detected               | p-Isopropyltoluene                  | Not Detected               |
| Butylbenzene (N- / tert- / sec-)                                              | Not Detected               | N-propylbenzene                     | Not Detected               |
| Carbon tetrachloride                                                          | Not Detected               | Naphthalene                         | Not Detected               |
| Chlorobenzene                                                                 | Not Detected               | Styrene                             | Not Detected               |
| Chloroform                                                                    | Not Detected               | 1,1,1,2-Tetrachloroethane           | Not Detected               |
| 2,4-Chlorotoluene                                                             | Not Detected               | Tetrachloroethene                   | Not Detected               |
| Dibromochloromethane                                                          | Not Detected               | Trichloroethene                     | Not Detected               |
| 1,2-Dibromo-3-chloropropane                                                   | Not Detected               | Trichloroethane (1,1,1- / 1,1,2-)   | Not Detected               |
| Dichlorobenzene (1,2- /1,3- /1,4-)                                            | Not Detected               | 1,2,3-Trichloropropane              | Not Detected               |
| Dichloroethane (1,1- / 1,2-)                                                  | Not Detected               | Trimethylbenzenes (1,2,4- /1,2,5-)  | Not Detected               |
| Dichloroethene (1,1- / cis-1,2- / trans-1,2-)                                 | Not Detected               | Toluene                             | Not Detected               |
| Dichloropropane (1,2- / 1,3-)                                                 | Not Detected               | Trichlorobenzenes (1,2,4- / 1,2,3-) | Not Detected               |
| 1,1-Dichloropropene                                                           | Not Detected               | o-Xylene                            | Not Detected               |

| Compounds                                                                                                                                                                                                                               | Whole device (µg/g or ppm) | Compounds                        | Whole device (µg/g or ppm) |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|----------------------------------|----------------------------|
| Ethanol, 2-trimethylsilyl                                                                                                                                                                                                               | 2.91 max.                  | m- & p-Xylene                    | Not Detected               |
| Total volatiles                                                                                                                                                                                                                         | 2.91                       |                                  |                            |
| Extractable*2 - Chemicals that can be released by breast implants following soaking in water                                                                                                                                            |                            |                                  |                            |
| Cyclic siloxanes (D4, D5, D6...)                                                                                                                                                                                                        | Not Detected               | Linear siloxanes (L3, L4, L5...) | Not Detected               |
| Total extractables                                                                                                                                                                                                                      | < 24.2                     |                                  |                            |
| *Not Detected means that the level of the individual volatile was below the quantitation limit of the test method. The quantitation limit was 0.86 µg/g for m- & p-Xylene and 2,4-Chlorotoluene, and 0.44 µg/g for all other compounds. |                            |                                  |                            |
| *Not Detected means that the level of the individual extractable was below 24.2 µg/g, the quantitation limit of the test method.                                                                                                        |                            |                                  |                            |

Table 3. Heavy Metals Found in Breast Implants

| Heavy metals                                                                                                               | Concentration (ppm) | Heavy metals | Concentration (ppm) |
|----------------------------------------------------------------------------------------------------------------------------|---------------------|--------------|---------------------|
| Arsenic                                                                                                                    | Not Detected        | Chromium     | Not Detected        |
| Cadmium                                                                                                                    | Not Detected        | Copper       | Not Detected        |
| Lead                                                                                                                       | Not Detected        | Tin          | Not Detected        |
| Mercury                                                                                                                    | 10 max.             | Nickel       | Not Detected        |
| Vanadium                                                                                                                   | Not Detected        | Platinum     | Not Detected        |
| Molybdenum                                                                                                                 | Not Detected        | Zinc         | 0.25 max.           |
| Selenium                                                                                                                   | Not Detected        | Magnesium    | 0.41 max.           |
| Cobalt                                                                                                                     | Not Detected        | Manganese    | Not Detected        |
| Antimony                                                                                                                   | Not Detected        | Thallium     | Not Detected        |
| Barium                                                                                                                     | Not Detected        |              |                     |
| Not Detected means that the level of the individual element was below 0.25 ppm, the quantitation limit of the test method. |                     |              |                     |

3. Warnings, Precautions to be Taken

Mammography

Women who have breast implants should continue to talk to their physician to carry out normal checks to detect breast cancer. The implant can hinder the detection of early breast cancer through mammography by concealing some underlying breast tissue. The implant can also compress overlying tissue which can 'hide' suspicious lesions in the breast. Please make sure you share all information you know about your breast implants with your healthcare provider so that they can adapt the mammography pressure accordingly. You should request radiologists who are experienced with the most current techniques and equipment for imaging breasts with implants. You should inform them of the presence, type and placement of implants and the orientation marks on the anatomical implant. The presence of a breast implant may make screening mammography more difficult and may lead to additional mammography views.

Detection of Breast Cancer

The presence of breast implants could delay the detection of breast cancer by self-exam. For this reason, patients must be informed that they need to consult a surgeon or doctor for appropriate medical monitoring. Patients should also get regular breast cancer screening.

Possible Interferences of the Implant in MR Environment / Security Scanners.

Please note that Eurosilicone's Gel-filled Breast Implants are all made from medical implant grade silicone materials which can be used with MRI Scan. No precaution to be taken.

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GC Aesthetics®

A Confident Choice for LifeTM  
Global Consolidated Aesthetics Limited, registered in Ireland with number 450181.  
Registered office:  
Suite 601, Q House, Furze Road, Sandyford Industrial Estate, Dublin 18, Ireland.  
© 2022 GC Aesthetics.® 904651.

GC Aesthetics® is a registered trademark  
of Global Consolidated Aesthetics Limited.

For more information about our products please visit:  
gcaesthetics.com or  
email: info@gcaesthetics.com